

PRESCOTT'S, INC.

Surgical Microscope Specialists

OmniPlus/Featherlite

Technical Data

| | |
|----------------------|--|
| Magnification | <i>5 Step Apochromatic magnification changer (0.4, 0.6, 1.0, 1.6, 2.5)</i> |
| Focus | <i>10 mm Fine Focus Objective Lens: 200mm, 250mm, 300mm, 400mm</i> |
| Illumination | <i>60 watt LED 12v, 150w Halogen</i> |
| Binocular | <i>F170 straight binocular F170 0-180 inclinable binocular</i> |
| Balancing | <i>Mechanical</i> |
| Asepsis | <i>Drapes, sterilizable caps</i> |
| Video | <i>Optional 4K recording and still image capture</i> |



ENGINEERED FOR MOBILITY

OmniPlus/Featherlite

Single Custom Pelican Case for Easy Transport



LED Illumination

- Superior brightness approximating natural light, CRI 92
- Long lasting bulb, estimated life 50,000 hours
- 5700K Color Temperature
- Turret Design on LED will accept multiple cable sizes
- Silicon based liquid fiber optic cable for improved light transmission



Optional 4k Video

- Integrated 4K Streaming video and still image capturing
- 64G memory
- Integrated camera can be paired with phone, tablet, or remote control
- Designated HDMI cable for monitor display
- USB 3.0/SD Card Reader

800-438-3937
www.surgicalmicroscopes.com

PRESCOTT'S, INC.
Surgical Microscope Specialists