



OMNIPLUS OSM200


DENTISTRY

**SURGICAL & CLINICAL
DOCUMENT IN 4K**



**NATIONWIDE COVERAGE
SEE BETTER**

PRESCOTT'S, INC.
Surgical Microscope Specialists



OMNIPLUS OSM200

For better results, every day.

Incorporate everything needed to fully concentrate on the procedure at hand.

TECHNICAL DATA

Magnification	5-step apochromatic magnification changer
Focus	Objective lens 200-350mm options Fine focus or varioscope
Main Lens	f=200-350mm options Apochromatic
Binocular	0-180 inclinable with ergonomic angled wedge; 10x or 12.5x eyepieces
Illumination	60w LED (5500K) paired with liquid light guide; optimal transmission; orange filter
Video	4K documentation, still and live image
Suspension	Rolling stand, wall and ceiling mount
Asepsis	Drapes, drape ring adapter, sterilizable caps
Accessories	Extension arm, rotational plate, MORA, onboard monitor


SAVE ON SPACE

With rolling stand, wall mount, and ceiling mount options.



DENTAL APPLICATIONS

- Endodontics
- Restorative
- Aesthetic
- Periodontics
- Implantology
- Educational use



OPTIMAL CLARITY SUPERIOR FOCUS & MAGNIFICATION

- 0-180 inclinable binoculars
- Ergonomic wedge
- 10x and 12.5x eyepieces
- High precision optics with broadband anti-reflexive coating
- Apochromatic 5-step magnification changer
- Lens options: 250-350mm focusable and varioscope



SEE BETTER LED ILLUMINATION & LLG

- Superior light output allows for better viewing of deep and narrow canals
- Liquid light guide (LLG) guarantees maximum light transmission
- 60w LED ideal for UV curing
- Color Rendering Index rating of 92
- Resembles natural daylight at 5700K color temperature
- 50,000 hours lifespan
- Easily adjust light intensity
- Accepts multiple cable sizes



DOCUMENT IN 4K STILL & LIVE IMAGE

- Integrated 4k camera
- Education and lecture uses
- SD memory card
- Wireless remote option




CUSTOMIZABLE TO MEET YOUR NEEDS

- Wall mount, ceiling mount, rolling stand options
- Long or short swing arm
- Integrated or non-integrated video
- Onboard monitor with mount
- Extension arm
- MORA attachments
- F170/F260 binoculars

NATIONWIDE COVERAGE

PRESCOTT'S, INC.
Surgical Microscope Specialists



CONTACT PRESCOTT'S TODAY TO GET IN TOUCH WITH YOUR EXCLUSIVE REPRESENTATIVE!

ASK ABOUT TRADE-IN CREDIT FOR YOUR OLD INVENTORY!

“ The Prescott's OmniPlus/OSM200 is an excellent office microscope for the ENT surgeon. It is very affordably priced and functions exceptionally well. I give Prescott's my highest recommendation.”

— Dr. Kuperan, *Houston Physicians Hospital*

“ I like the ease of use. Everything is built in and I can record very easily.”

— Yaara Berdan, *Herman Ostrow School of Dentistry of USC*

“ Scope function is easy and simple. Visual acuity is top notch, and LED lighting is exactly what we needed.”

— Dr. Montano, *McAllen Surgical Specialty Center*

Prescott's, Inc.
18940 Microscope Way
Monument, CO 80132

Phone: (800) 438-3937
prescott@surgicalmicroscopes.com
www.surgicalmicroscopes.com



EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands