



HD equipment that goes where it's needed.

Sony's HD on the Go is the mobile solution for multiple ORs.

Now there's an HD video solution for your operating rooms that is both flexible and economical. Sony's HD on the Go is a mobile cart that can be moved to wherever it's needed, stocked with all the HD video equipment necessary for today's ORs, and even customized to your specifications. When budgets are an issue, HD on the Go offers a simple way to equip your ORs without purchasing multiple sets of equipment.



Benefits:

- Plug & Play sharing between multiple ORs
- No pre-wiring needed
- Flexible configurations for multiple applications
- Compatible with all camera systems
- Built-in isolation transformer available

HD equipment that does multiple duty.

HD on the Go will transform your OR experience. It's flexible and can easily be fully outfitted with the latest Sony HD video equipment customized to your needs. And the whole set-up is lightweight, portable, and easy to position in any room.

LMD-1530MD 15" * HD Medical Grade Widescreen Monitor. Equipped with a 10-bit signal processor to reproduce lifelike images with smooth gradation and accurate color consistency.

Medical Grade Keyboard. Features a completely sealed water-resistant robust rubber casing which can be washed down to prevent the spread of bacteria.

UP-DR80MD Medical Grade Printer. Full-page digital dye-sublimation color printer with outstanding image quality and accurate color reproduction.

HVO-1000MD Medical Blu-ray Recorder. Simultaneous recording on a 320GB internal hard drive and various external media eliminates the need for time-consuming secondary image duplication. Also has the ability to support live image streaming.

LMD-2451MD 24" * HD Medical Grade Widescreen Monitor. Offers a superb quality WUXGA (1920 x 1200) LCD panel, 10-bit signal processing, Sony ChromaTRU™ color balance.

PMW-10MD HD Video Camera. Equipped with 3 1/2-inch (1920 x 1080) Exmor™ CMOS sensors, each delivering over 2 million pixels.

The configuration shown is only one of the many ways HD on the Go can be customized to your needs.

When your budget is limited but your need for HD is not, use HD on the Go to expand the HD capacity in your ORs. It's the efficient and affordable way to put superb quality HD right where it's needed at a fraction of the cost of multiple systems.

Visit sony.com/hdforsurgery for more information.

*Viewable area, measured diagonally.

Sony Electronics Inc.
1 Sony Drive
Park Ridge, NJ 07656
sony.com/hdforsurgery

©2011 Sony Electronics Inc. All rights reserved. Reproduction in whole or in part without written permission is prohibited. Features and specifications subject to change without notice. Sony, make.believe, Blu-Ray, Chroma-Tru, and Exmor and their respective logos are trademarks of Sony. All other trademarks are properties of their respective owners.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician or other appropriately licensed medical professional. CAUTION: See product labeling for indications, contraindications, warnings, cautions, and directions for use.